



## Strong correlation between in vitro microfluidic skin model comprising precise sebum and sweat excretion with in vivo studies for makeup lasting claim substantiation

Débora Freitas<sup>1</sup>, Pascale Barlier<sup>2</sup>, Kaeshna Rughoobur<sup>3</sup>, Muaawiyah Adaam Domun<sup>3</sup>, Dylan Jodun<sup>3</sup>, Fabrice Monti<sup>1\*</sup>

- <sup>1</sup> Microfactory SAS, Paris, France
- <sup>2</sup> Cos&Co SAS, Versailles, France
- <sup>3</sup> Centre International de Développement Pharmaceutique (CIDP), Phoenix, Mauritius
- \* Corresponding author: Fabrice Monti, +33662720461, fabrice.monti@microfactory.eu

#### Introduction

Long lasting effect and mattifying effect have become basics requirements for foundation products. In vivo studies are the gold standard method for testing makeup products, by applying them to

multiple volunteers. Different methods can be chosen depending on the desired claims. However, this process takes several months to discriminate between makeup formulations. To find an alternative to volunteers, we have developed a fast and innovative in vitro method.

The Microfluidic Universal Skin (MUS) is an innovative microfluidic skin model that mimics the mechanism of sweat and sebum excretion. MUS is therefore ideal for evaluating the in vitro efficacy of makeup products regarding their resistance to sebum and sweat excretion.

To assess the effectiveness of MUS, we compared the results of in vivo tests conducted on volunteers with those obtained in vitro with MUS.

#### **Materials & Methods**

In vitro protocol	In vivo protocol
<ul> <li><u>Controlled parameters:</u></li> <li>Artificial sweat, pH=5.5 [1] &amp; artificial sebum, viscosity = 56 mPa.s [2]</li> <li>Microfluidics system designed to mimic <i>in vitro</i> sebum and sweat excretion under controlled parameters such as flow injection</li> <li>Replication of human skin characteristics: gland density, roughness and hydrophobicity</li> </ul>	<ul> <li>Inclusion criteria:</li> <li>12 women selected</li> <li>Combination skin</li> <li>From 18 to 65 years old</li> <li>Phototype: II - IV</li> </ul>
<ul> <li>3 liquid foundations tested</li> <li>3,8 mg/cm<sup>2</sup></li> </ul>	<ul> <li>3 liquid foundations tested</li> <li>2 mg/cm<sup>2</sup></li> </ul>
<ul> <li>Methodology:</li> <li>Standard images, parallel polarized images and crossed polarized images using MUS imaging system</li> <li>Continuous analysis (1 image every 2 minutes)</li> </ul>	<ul> <li><u>Methodology:</u></li> <li>4 timepoints (T<sub>0</sub>, T<sub>immediate</sub>, T<sub>12h</sub>, T<sub>24h</sub>)</li> <li>Clinical scoring by trained experimenter</li> <li>Front view images using the VISIA-CR® (standard, parallel polarized and crossed polarized images)</li> <li>Image analysis by Newtone</li> </ul>

#### lage analysis by Newtone

#### **Results & Discussion**

#### 1) Color variation

*In vivo*: The clinical scoring study shows a similar trend between 12 and 24 hours. However, due to limitations in the human eyes gradation capabilities, the image



analysis study is more accurate thanks to a multi-parametric analysis. In fact, the color variation of the products in image analysis shows the same trend at 12 and 24 hours. The product C has a higher color variation compared to products A and B. The product A is evaluated as the best product of all three in terms of longlasting effect.





The MUS model is an accelerated sweat and sebum excretion process. Therefore, 1 hours of in vitro MUS corresponds to

12 hours in vivo, and 2 hours in vitro corresponds to 24 hours in vivo. The three products are impacted differently by the two fluids, sebum and sweat. The cumulated color variation due to sebum and sweat is studied to predict their lasting effect over 12 hours and 24 hours.

Sebum influence Sweat influence





The gloss parameter is used to determine resistance to sebum and sweat and to assess if a product has a mattifying effect. The gloss variation of the products on *in vivo* image analysis shows that product B is less resistant than the other two products C and A. MUS is able to predict the same trends as in vivo analyses. Product A has a low gloss variation, in contrast products C and B, which have a higher gloss variation.



**Discussion:** Products A, B and C could be evaluated with two methods, *in vivo* analysis and the *in vitro* MUS model. These two models predicted the same resistance to fluids whether in gloss or color variation. By these two analyses, product A is the most resistant. The MUS model permit to predict the performance of products but also to provide informations on the impact of sebum and sweat.

## Conclusions





# References

[1] C.J. Harvey, R.F. Lebouf, A.B. Stefaniak. Formulation and stability of a novel artificial human sweat under conditions of storage and use (2010) [2] E.O. Butcher and A. Coonin. *The physical properties of human sebum* (1949)

In this study, we investigated the long lasting effect (color evolution & gloss evolution) of three liquid foundations using three distinct methodologies:

- 1. The clinical scoring method represents the subjective perception of a trained human eye and is extensively utilized by professionals as it provides insights into consumer perception.
- 2. Image analysis involves sophisticated algorithms with enhanced discriminatory power.
- 3. MUS replaces human volunteers with a skin model capable of mimicking sweat and sebum excretion. This skin model is combined with sophisticated image analysis to assess color and gloss variations.

The results obtained from these three techniques were in excellent agreement. To conclude, MUS is an *in vitro* model correlated to *in vivo* for the evaluation of makeup products.





©2023 Copyright by Author for the Sociedad Española de Químicos Cosméticos (SEQC)'s 33rd IFSCC Congress. Unauthorized reproduction and/or distribution prohibited without prior written consent of copyright holder and the SEQC.

33<sup>rd</sup> IFSCC Congress, 4-7 September 2023, Barcelona.

#### **Rethinking Beauty Science**